



RELAX RANGE OF PRESSURE CUSHIONS WARRANTY AND INFORMATION

*RELAX wheelchair cushions are guaranteed against manufacturing defects for 2 years from date of manufacture.

The inner cushion should be inspected regularly – if damaged it should not be used.

The cover if damaged should not be used.

It is very important to reposition the user on a regular basis.

To ensure maximum shelf life, it is recommended to store the cushion flat in its box.

If your RELAX wheelchair cushion has been subjected to extreme temperatures (-20°C/-4°F or +40°C/104°F) please allow to revert to room temperature before use

*Pressure mapping information available on request.

*Each RELAX cushion is recommended to be used with the cover provided.

When fitting the cushion cover ensure the back of the cushion is positioned to the back of the cover.

It is important to clean the cushion and cover between users to avoid cross contamination.

Cleaning instructions

Outer Cover: Machine Wash 60°C (140°F). Dry Flat. Do Not Tumble Dry.

Inner Cover: The Interior Cover can be cleaned using a damp cloth with a mild disinfectant and should be cleaned regularly.

Disposal instructions

Any disposal practice must be in compliance with all local and national laws and regulations.

All packaging materials used for your RELAX cushion can be recycled.

If returning the product, the cover should be washed and the inner cushion disinfected before being placed in a sealed plastic bag.